

# Storage options

## MAXBIN™ STORAGE SYSTEM

The MaxBin system provides high-capacity storage for medication and supplies, in your choice of electronic locking or non-locking bins.

- ▶ Higher capacity, more flexible storage system
- ▶ 60-74% more capacity per row than standard bins
- ▶ Full width drawers with soft-close slides
- ▶ Dual locking option for controlled medications and witness authentication
- ▶ Remote security management

## STANDARD BIN SYSTEM

Trio is configurable with a selection of standard drawers and bins in electronic locking or non-locking. Mix and match to meet your facility's workflow, storage and security requirements. Bins are available in gray or blue.

- ▶ Electric locking or non-locking drawers and bins
- ▶ Available in 2, 4 or 6 tier storage options
- ▶ Arrange your storage with a selection of bin and drawer configurations
- ▶ Available colors include NEW gray or blue
- ▶ Remote security management




*Trio with  
Standard Bin System*

# Standard Bin System

## Standard Bin Configurations

Configure the Trio workstation with standard drawer and bin options in your choice of electronic locking or non-locking to meet specific workflow, storage and security needs. Available in gray or blue.

### Preconfigured Bin Modules



### Standard Bin Dimensions



**3.875" Small Bin**  
3.875" w x 2.375" h x 10.5" d  
(9.84 x 6.03 x 26.67 cm)



**6.375" Medium Bin**  
6.375" w x 2.375" h x 10.5" d  
(16.19 x 6.03 x 26.67 cm)



**12.75" Wide Bin**  
12.75" w x 2.375" h x 10.5" d  
(32.38 x 6.03 x 26.67 cm)



**6.375" Medium Deep Bin**  
6.375" w x 4.9375" h x 10.5" d  
(16.19 x 12.54 x 26.67 cm)



**12.75" Wide Deep Bin**  
12.75" w x 4.9375" h x 10.5" d  
(32.38 x 12.54 x 26.67 cm)



*Envelope Drawer available in blue only.*

**Trio with Standard Bins**

# MaxBin™ System

## MaxBin High Capacity Configurations

The MaxBin system provides high-capacity storage for medication and supplies, in your choice of electronic locking or non-locking bins. The MaxBin's unique configurable design provides options for simple 1-tier up to expansive 5-tier designs.



**Trio with MaxBin System**

### Sample Tier Configurations



### MaxBin Dimensions



**4" Small Bin**  
4" w x 2.66" h x 11.86" d  
(10.16 x 6.75 x 30.12 cm)



**5.5" Medium Bin**  
5.5" w x 2.66" h x 11.86" d  
(13.97 x 6.75 x 30.12 cm)



**8.6" Wide Bin**  
8.6" w x 2.66" h x 11.86" d  
(21.59 x 6.75 x 30.12 cm)



**16.6" Extra Wide Drawer**  
16.6" w x 2.66" h x 11.86" d  
(42.16 x 6.75 x 30.12 cm)



**4" Small Deep Bin**  
4" w x 6.15" h x 11.86" d  
(10.16 x 15.62 x 30.12 cm)



**5.5" Medium Deep Bin**  
5.5" w x 6.15" h x 11.86" d  
(13.97 x 15.62 x 30.12 cm)



**8.6" Wide Deep Bin**  
8.6" w x 6.15" h x 11.86" d  
(21.59 x 15.62 x 30.12 cm)



**16.6" Extra Wide Deep Drawer**  
16.6" w x 6.15" h x 11.86" d  
(42.16 x 15.62 x 30.12 cm)