

HiCare E

Combines advanced ergonomics, durable construction, and premium mobility for point-of-care computing.



Clinical Features

- Adjustable LED keyboard, storage, and task lighting
- Antimicrobial touch surfaces
- 8' hospital-grade power cord
- Easy-to-clean, disinfectant-compatible surfaces
- Raised-edge worksurface with Lexan cover
- Multiple accessory mounting locations
- User-friendly controls

Ergonomic Features

- 270° wraparound handles
- Patented ergonomic adjustable push handles
- One-touch worksurface height adjustment
- One-touch LCD height adjustment
- Patented on-demand electronic steering assist
- Rubberized footrest on each leg
- Slide-out keyboard tray with wrist rest and auto-leveling tilt

Technical Features

- AC power distribution w/6 AC receptacles
- Cart management software
- CE and UL certified
- Easily accessible keyboard and mouse connections
- Electronic locking PC compartment
- Thermally controlled PC compartment cooling fan
- Tool-less battery replacement
- Ships fully assembled
- 7-port USB hub
- External USB charging port

General Specifications

- Base footprint: 19" W × 20" D
- Overall footprint: 21.6" W × 23.6" D
- Usable work surface: 17" W × 14.3" D
- PC compartment: 17.1" W × 3.2" H × 12.1" D
- Keyboard tray: 18.2" W × 1.6" H × 7.5" D
- Keyboard height adjustment: 26.9" to 42.0"
- Monitor height: 42.0" to 62.4"
- Monitor adjustment range: 5.3"
- Push handle height: 30.3" to 45.4"
- Starting weight: 112 lbs.

- Push/pull force: 1.8 lbs.
- Maximum LCD weight: 20 lbs.
- 3-year standard warranty

Height Adjustment

- No storage: 32.5" - 47.6"
- With 2-tier storage: 32.5" - 47.6"
- With 4-tier storage: 34.1" - 47.6"
- With 6-tier storage: 40.7" - 47.6"

Power System

- 8' coiled dual hook hospital-grade power cord
- 150-watt inverter with 20A charger
- 120VAC true sine wave output
- 50Ah lithium iron phosphate (LiFePO₄) battery
- Charge time from empty: 2.75 hours
- Battery life expectancy: 5+ years