

Howard HiCare Cart Comparison

For upgrades and add-ons, contact your
Convergent representative.

	HiCare Lite	HiCare E
General Specs		
Base	19"W x 20"D	
Worksurface	17"W x 14.3"D	
PC Compartment	18.2"W x 12.5"D x 3.2"D	
Battery	50ah lithium iron phosphate	
Height Adjustment		
One-touch powered work surface	32.5" to 47.6"	
One-touch powered LCD arm		5.25" Range
Fixed height LCD arm	o	
Clinical Specifications		
Raised edge work surface with document cover	o	
5" Dual wheel casters with 2 locking casters	o	
Antimicrobial touch surfaces	o	
Multi-brightness level LED for keyboard, work surface, and storage	o	
Ergonomic Specifications		
Patented ergonomically adjustable push-handles	o	
Rubberized footrest	o	
-/+ Tilting slide-out keyboard tray with integrated wrist rest	o	
VESA Mount with tilt, swivel, and rotate	o	
Central control panel	o	
270 Degree wrap-around handle	o	
Patented on-demand steer assist technology	o	
Technical Specifications		
PC Compartment cooling fan	o	
AC Power distribution with 6 AC receptables and 7-port USB hub	o	
Keyboard and mouse connections	o	
Detent & key locking storage (up to 2 rows of drawers)	o	
Key locking PC compartment	o	
6' Medical-grade power cord	o	
Electronic locking storage with independent key bypass (up to 6 rows of drawers)		o
Electronic locking PC compartment		o
8' Medical-grade power cord		o
External USB charging port		o
Ships fully assembled	o	
UL & CE Certified	o	



Trusted partner. Trusted solutions.

innovative-medical.com
114info@convergint.com
800.968.3135

convergint® | innovative
medical