

# Howard HiCare Cart Comparison

For upgrades and add-ons, contact your  
Convergent representative.

|                                                                                  | HiCare Lite             | HiCare E    |
|----------------------------------------------------------------------------------|-------------------------|-------------|
| General Specs                                                                    |                         |             |
| Base                                                                             | 19"W x 20"D             |             |
| Worksurface                                                                      | 17"W x 14.3"D           |             |
| PC Compartment                                                                   | 18.2"W x 12.5"D x 3.2"D |             |
| Height Adjustment                                                                |                         |             |
| One-touch powered work surface                                                   | 32.5" to 47.6"          |             |
| One-touch powered LCD arm                                                        |                         | 5.25" Range |
| Fixed height LCD arm                                                             | o                       |             |
| Power Options                                                                    |                         |             |
| 38.4 Ah Battery, 491 whr                                                         | o                       |             |
| 51.2 Ah Battery, 655 whr                                                         |                         | o           |
| Clinical Specifications                                                          |                         |             |
| Raised edge work surface with document cover                                     | o                       |             |
| 5" Dual wheel casters with 2 locking casters                                     | o                       |             |
| Antimicrobial touch surfaces                                                     | o                       |             |
| Multi-brightness level LED for keyboard, work surface, and storage               | o                       |             |
| Ergonomic Specifications                                                         |                         |             |
| Patented ergonomically adjustable push-handles                                   | o                       |             |
| Rubberized footrest                                                              | o                       |             |
| -/+ Tilting slide-out keyboard tray with integrated wrist rest                   | o                       |             |
| VESA Mount with tilt, swivel, and rotate                                         | o                       |             |
| Central control panel                                                            | o                       |             |
| 270 Degree wrap-around handle                                                    | o                       |             |
| Patented on-demand steer assist technology                                       | o                       |             |
| Technical Specifications                                                         |                         |             |
| PC Compartment cooling fan                                                       | o                       |             |
| AC Power distribution with 6 AC receptables and 7-port USB hub                   | o                       |             |
| Keyboard and mouse connections                                                   | o                       |             |
| Detent & key locking storage (up to 2 rows of drawers)                           | o                       |             |
| Electronic locking storage with independent key bypass (up to 6 rows of drawers) |                         | o           |
| Electronic locking PC compartment                                                |                         | o           |
| Key locking PC compartment                                                       | o                       |             |
| 8' Medical-grade power cord                                                      |                         | o           |
| 6' Medical-grade power cord                                                      | o                       |             |
| External USB charging port                                                       |                         | o           |
| Ships fully assembled                                                            | o                       |             |
| UL & CE Certified                                                                | o                       |             |



Trusted partner. Trusted solutions.

innovative-medical.com  
 114info@convergint.com  
 800.968.3135

convergint® | innovative  
 medical