



HI-CARE LITE LCD

FEATURES AND SPECIFICATIONS



CLINICAL FEATURES

- Raise-Edge Work Surface w/ Lexan Cover
- 5" Dual Wheel Caster w/Two Locking Casters
- Anti-Microbial Agents in Touch Surfaces
- Multi Bright Level LED Keyboard Lighting
- Multi Bright Level Red & White LED Work Surface Light
- LED Storage Lighting

ERGONOMIC FEATURES

- Rubberized foot rests on each side of the cart allow for a more relaxed stance and relieves fatigue
- Slide-out keyboard tray with integrated wrist rest is auto-leveling with a 15-degree negative-tilt and a positive-tilt up to 6-degrees
- One-touch worksurface height adjustment allows for a full sit/stand range of motion
- 270° VESA mount to allow full view swivel when maneuvering cart
- Ergonomically positioned optional storage module
- Open architecture at wheel base allows end user to maintain their normal walking stride
- 270° Warp around handles allows for enhanced ergonomics
- Equipped with patented on-demand steer assist technology (SAT)
- Fixed LCD Arm Height Adjustment
- Patented Ergonomic Adjustable Push Handles

TECHNICAL FEATURES

- Ships Fully Assembled
- Tool-Less Battery Access and Replacement
- Easily Replaceable Power Module
- AC Power Distribution w/6 AC Receptacles
- Easily Accessible Keyboard and Mouse Connections
- Thermally Controlled PC Compartment Cooling Fan
- UL and CE Certified
- Storage Options Available
- 6' Hospital Grade Power Cord
- Keyed Locking PC Compartment

SPECIFICATIONS

GENERAL

- Base: 19" W x 20" D
- Work Surface: 17" W x 14.3" D
- PC Compartment: 17.1" W x 3.2" H x 12.5" D
- Keyboard Tray: 18.2" W x 1.6" H x 7.5" D
- Keyboard Tilt: Negative 5° / Positive 6°
- Max Monitor Weight: 20 lbs.

HEIGHT ADJUSTMENTS

- Work Surface: 32.5" - 47.6" (no storage/2 Tier storage)
- Keyboard: 26.9" - 42.0"
- Handles: 30.3" - 45.4"
- Monitor: 42" - 62.4"

POWER SYSTEM

POWER MAX 150 POWER MODULE

- 150 W
- 20A Charger
- 120VAC True Sine Output

LITHIUM IRON PHOSPHATE BATTERY

- 38.4AH, 491Whr

STANDARD WARRANTY

CART

- 3-Year Electrical
- 3-Year Mechanical

BATTERY

- 3-Year Defect
- 3-Year Performance

Distributed by:



convergent[®] | innovative medical

www.innovative-medical.com | 800.968.3135 | 114info@convergent.com



HI-CARE LITE

STORAGE CONFIGURATION OPTIONS



2 TIER DETENT OPTIONS



(1x) Two-Tier Drawer



(2x) One-Tier Drawers



(1x) One-Tier Drawer
(1x) Two-Tier Hang Down Drawer

2 TIER KEYED OPTIONS



(1x) Two-Tier Drawer



(2x) One-Tier Drawers



(1x) One-Tier Drawer
(1x) Two-Tier Hang Down Drawer

2 TIER KEYPAD LOCKING OPTIONS



(1x) Two-Tier Drawer



(2x) One-Tier Drawers



(1x) One-Tier Drawer
(1x) Two-Tier Hang Down Drawer

LOCKING OPTIONS

Locks

- Detent
- Keypad
 - USB Powered or AAA Batteries
 - Up to 28 4-Digit Codes
- Keypad
 - Includes 2 Keys Per Cart

STORAGE OPTIONS

STORAGE

- 2 Tiers of Configurable Storage
- 3 Storage Configuration Options
- Storage Drawers (Dividers Included)

DRAWER DIMENSIONS



1 Tier Drawer: 12.2" W x 2.5" H x 10.9" D



2 Tier Drawer: 12.2" W x 5.4" H x 10.9" D



2 Tier Hang Down Drawer: 12.2" W x 5.4" H x 10.9" D



HI-CARE LITE LCD ACCESSORIES OPTIONS



24" Monitor



Dual Monitor Arm



Wireless Scanner Mount



Tether Scanner Mount



External Antenna Mount



Quick Release



Card Reader Bracket¹



Cup Holders



Zebra Printer Bracket¹



Washable Keyboard



Washable Mouse



RFID Reader Bracket



Integration Kit



7-Port USB Hub



XL Worksurface



1L Hand Sanitizer Bracket



2.2qt Sharps Bracket



Glove Box Bracket



Hand Sanitizer Bracket



Mini-Waste



Sani Wipe Bracket



E Signature Pad



Handle w/ Basket



XL Device Shelf



*Accessory options may vary between carts

*Mounting location are not available for all Howard carts

¹Secondary mounting location available