

Clinical Powered Carts

HI-CARE COMPARISON



	HI-Care Lite LCD	HI-Care Lite Laptop	HI-Care E
General Specs			
Base	19"W x 20"D		
Worksurface	17"W x 14.3"D		
PC Compartment	18.2"W x 12.5"D x 3.2"H		
Height Adjustment			
One-Touch Powered Work Surface (w/ No Storage)	32.5" to 47.6"		
One-Touch Powered LCD Arm (Electronic)			5.25" Range (Electronic)
LCD Arm (Fixed height)	✓		
Power Options			
38.4 Ah Battery, 491 Whr	✓	✓	
51.2 Ah Battery, 655 Whr			✓
Clinical Specifications			
Raised Edge Work Surface w/ Lexan Cover	✓	✓	✓
5" Dual Wheel Casters w/ Two Locking Casters	✓	✓	✓
Medical-Grade Touch Surfaces	✓	✓	✓
Multi-Brightness Level LED Keyboard Lighting	✓	✓	✓
Multi-Brightness Level Red & White LED Worksurface Light	✓	✓	✓
LED Storage Lighting	✓	✓	✓
Ergonomic Specifications			
Patented Ergonomically Adjustable Push-Handles w/ Soft-Touch Grip	✓	✓	✓
Rubberized Footrest	✓	✓	✓
-/+Tilting Slide-Out Keyboard Tray w/ Integrated Wrist Rest	✓	✓	✓
One-Touch Powered Work Surface Height Adjustment	✓	✓	✓
Ergonomically Positioned Storage Module	✓	✓	✓
VESA Mount has Tilt, Swivel, and Rotation Adjustments	✓	✓	✓
Centralized Control Panel	✓	✓	✓
270 Degree Wrap-Around Handle	✓	✓	✓
Patented On-Demand Steering Assist Technology	✓	✓	✓
Technical Specifications			
Shipped Fully Assembled	✓	✓	✓
Easily Replaceable Power Module and Battery	✓	✓	✓
Thermally Controlled Variable Speed PC Compartment Cooling Fan	✓	✓	✓
AC Power Distribution w/ 6 AC Receptacles	✓	✓	✓
Easily Accessible Keyboard and Mouse Connections	✓	✓	✓
Detent Only Storage	✓ (see storage guide)	✓ (see storage guide)	
Keyed Locking Storage	✓ (see storage guide)	✓ (see storage guide)	
Electronic Locking Storage w/ Independent Key Bypass			✓ (see storage guide)
Electronic Locking PC Compartment			✓
Keyed Locking PC Compartment	✓	✓	
8' Hospital-Grade Power Cord			✓
6' Hospital-Grade Power Cord	✓	✓	
External USB Charging Port			✓
UL and CE Certified	✓	✓	✓

