

HiCare Comparison

	HiCare Lite	HiCare E
General Specs		
Base	19"W x 20"D	
Worksurface	17"W x 14.3"D	
PC Compartment	18.2"W x 12.5"D x 3.2"D	
Height Adjustment		
One-Touch Powered Worksurface	32.5" to 47.6"	
One-Touch Powered LCD Arm (Electronic)		5.25" Range
Fixed Height LCD Arm	x	
Power Options		
38.4 Ah Battery, 491 Whr	x	
51.2 Ah Battery, 655 Whr		x
Clinical Specifications		
Raised Edge Worksurface with Lexan Cover	x	x
5" Dual Wheel Casters with 2 Locking Casters	x	x
Medical-Grade Touch Surfaces	x	x
Multi-Brightness Level LED Keyboard, Worksurface, and Storage Lighting	x	x
Ergonomics Specifications		
Patented Ergonomically Adjustable Push-Handles with Soft-Touch Grip	x	x
Rubberized Footrest	x	x
-/+ Tilting Slide-Out Keyboard Tray with Integrated Wrist Rest	x	x
One-Touch Powered Worksurface Height Adjustment	x	x
VESA Mount with Tilt, Swivel, and Rotation	x	x
Centralized Control Panel	x	x
270 Degree Wrap-Around Handle	x	x
Patented On-Demand Steering Assist Technology	x	x
Technical Specifications		
Ships Fully Assembled	x	x
Easily Replaceable Power Module & Battery	x	x
Thermally Controlled Variable Speed PC Compartment Cooling Fan	x	x
AC Power Distribution with 6 AC Receptacles	x	x
Easily Accessible Keyboard & Mouse Connections	x	x
Detent & Key Locking Storage (Up to 2 Rows of Drawers)	x	
Electronic Locking Storage with Independent Key Bypass (Up to 6 Rows of Drawers)		x
Electronic Locking PC Compartment		x
Key Locking PC Compartment	x	
8' Hospital-Grade Power Cord		x
6' Hospital-Grade Power Cord	x	
External USB Charging Port		x
UL & CE Certified	x	x

