


HiCare E

2-Tier Storage Configuration Options




01



(1x) Two-Tier Patient Bin

02




(2x) One-Tier Patient Bins

03



(4x) 6" Patient Bins

04



(6x) 4" Patient Bins

05



(3x) Two-Tier 4" Patient Bins

06



(2x) Two-Tier 6" Patient Bins

Additional Storage Options are Available Upon Request

Bin Dimensions:



1-Tier Bin 12" W x 2.25" H x 11"D



2-Tier Bin: 12" W x 5.25" H x 11"D



4" Bin: 4" W x 2.25" H x 11"D



6" Bin: 6" W x 2.25" H x 11"D



2-Tier 4" Bin: 4" W x 5.25" H x 11"D



2-Tier 6" Bin: 6" W x 5.25" H x 11"D

Storage Features:

- Electronic Locking
- Basic or Advanced Electronics Storage Cabinet Options Available
- Support 2000 Unique Codes
- Audit Logging
- Supports Independent, Second Level Locking for Securing Narcotics
- Open Drawer Alert (Advance Only)
- On/Off Storage Light
- On-Screen Drawer Configuration and Status via Software
- Independent Key Bypass
- 2, 4, or 6 Tiers of Configurable Storage
- 100+ Storage Configurations
- Patient Bins with Dividers
- Storage Drawers with Dividers
- Double Lock Narc Drawers (Optional)
- Quick Swap Cassette Exchange
- Adjustable Auto Locking Timeout

4-Tier Storage Configuration Options



01



(2x) Two-Tier Patient Bins

02



(2x) One-Tier Patient Bins
(1x) Two-Tier Patient Bin

03



(4x) 6" Patient Bins
(1x) Two-Tier Patient Bin

04



(6x) 4" Patient Bins

05



(4x) One-Tier Patient Bins

06



(8x) 6" Patient Bins

07



(12x) 4" Patient Bins

08



(3x) Two-Tier 4" Patient Bins
(1x) Two-Tier Patient Bin

09



(2x) Two-Tier 4" Patient Bins
(1x) Two-Tier Patient Bin

*Additional Storage Options are
Available Upon Request*

Bin Dimensions:



1-Tier Bin: 12" W x 2.25" H x 11"D



2-Tier Bin: 12" W x 5.25" H x 11"D



4" Bin: 4" W x 2.25" H x 11"D



6" Bin: 6" W x 2.25" H x 11"D



2-Tier 4" Bin: 4" W x 5.25" H x 11"D





2-Tier 6" Bin: 6" W x 5.25" H x 11"D

Storage Features:

- Electronic Locking
- Basic or Advanced Electronics Storage Cabinet Options Available
- Support 2000 Unique Codes
- Audit Logging
- Supports Independent, Second Level Locking for Securing Narcotics
- Open Drawer Alert (Advance Only)
- On/Off Storage Light
- On-Screen Drawer Configuration and Status via Software
- Independent Key Bypass
- 2, 4, or 6 Tiers of Configurable Storage
- 100+ Storage Configurations
- Patient Bins with Dividers
- Storage Drawers with Dividers
- Double Lock Narc Drawers (Optional)
- Quick Swap Cassette Exchange
- Adjustable Auto Locking Timeout

6-Tier Storage Configuration Options



Additional Storage Options are Available Upon Request

Bin Dimensions:



1-Tier Bin: 12" W x 2.25" H x 11"D



2-Tier Bin: 12" W x 5.25" H x 11"D



4" Bin: 4" W x 2.25" H x 11"D



6" Bin: 6" W x 2.25" H x 11"D



2-Tier 4" Bin: 4" W x 5.25" H x 11"D



2-Tier 6" Bin: 6" W x 5.25" H x 11"D

Storage Features:

- Electronic Locking
- Basic or Advanced Electronics Storage Cabinet Options Available
- Support 2000 Unique Codes
- Audit Logging
- Supports Independent, Second Level Locking for Securing Narcotics
- Open Drawer Alert (Advance Only)
- On/Off Storage Light
- On-Screen Drawer Configuration and Status via Software
- Independent Key Bypass
- 2, 4, or 6 Tiers of Configurable Storage
- 100+ Storage Configurations
- Patient Bins with Dividers
- Storage Drawers with Dividers
- Double Lock Narc Drawers (Optional)
- Quick Swap Cassette Exchange
- Adjustable Auto Locking Timeout