


HiCare Lite

Storage Configuration Options

2-Tier Detent Options

01



(1x) Two-Tier Drawer

02



(2x) One-Tier Drawer

03



(1x) One-Tier Drawer
(1x) Two-Tier Hang Down Drawer

2-Tier Keyed Locking Options

04



(1x) Two-Tier Drawer

05



(2x) One-Tier Drawer

06



(1x) One-Tier Drawer
(1x) Two-Tier Hang Down Drawer

2-Tier Keypad Locking Options

07



(1x) Two-Tier Drawer

08



(2x) One-Tier Drawers

09



(1x) One-Tier Drawer
(1x) Two-Tier Hang Down Drawer



Additional Storage Options are Available Upon Request


Bin Dimensions:



1-Tier Bin 12" W x 2.25" H x 11"D



2-Tier Bin: 12" W x 5.25" H x 11"D



2 Tier Hang Down Drawer:
2.2" W x 5.4" H x 10.9" D

Storage Features:

- Locking Options:
 - Detent
 - Keyed Locking (2 Keys Per Cart Included)
 - Keypad Locking
 - USB Powered or AAA Batteries
 - Up to 28 4-Digit Codes