

## DCT-6 – Powered Cart

### Computer Carts for Healthcare

*Proudly Made in the USA*

#### Superior Engineering

*Leading-edge design, production,  
and assembly on all products.*

#### Custom Configuration

*Building functionality that you want  
and need down to the last detail.*

#### Valued Experience

*Overall product performance and  
support at a valued price.*

#### Innovative Design

*Building better solutions that meet  
your specific needs and requirements.*



### Key Features

- Height Adjustable, Vertical Travel for Sit to Stand
- Electric Push Button for Stable Height Adjustment
- Durable and Reliable Construction for Ease of Access, Maintenance and Cleaning
- Secure, Adjustable Locking Compartment for Computer and Cable Management
- Small Footprint with Open Front for Effortless Transportation
- Long Lasting Lithium Iron Phosphate Power System with A/C Outlets and Remote User Interface
- Fleetview Management Software
- Laminate Work Surface with Raised Edge
- Easy Integration of Accessories
- Ergonomically Placed Footrest Pads
- (4) 5" Locking Casters
- Contoured Steering Handle
- Custom Configurations Available

**The Best Quality.  
The Best Value.**

Our reputation is for superior engineering, short turn-around times, strong customer service and a value proposition to retain your business for years to come.



## DCT-6 – Powered Cart

### Part Numbers:

**DCT-6LiA:** Secure Laptop Enclosure Cart with Powered Lift

**DCT-6LiA-LCD:** LCD Secure CPU Enclosure Cart with Powered Lift

**DCT-6LiA-P16V:** LCD Secure CPU Enclosure Cart with Extruded Pole VESA Mount

**DCT-6LiA-P16VA:** LCD Secure CPU Enclosure Cart with Extruded Pole VESA Articulating Arm Mount

**DCT-6LiA-P16DMVA:** LCD Secure CPU Enclosure Cart with Extruded Pole and Dual VESA Articulating Arm Mount



DCT-6LiA-P16V



DCT-6LiA-P16VA



DCT-6LiA-P16DMVA

### Technical Specifications:

- Work Surface Height:
  - Range: 14.00"
  - Standard Minimum: 33.5"
  - Standard Maximum: 47.5"
- Work Surface Area:
  - Standard: 19.75" x 18.0"
- Standard Wheels:
  - 5" High Performance Single Wheel Casters
- Base Footprint Dimensions:
  - 18.5" x 18.5"
- Locking PC Enclosure Dimensions:
  - 17.5" x 15.0" x 2.975"
- Keyboard Tray:
  - 18.0" x 9.0"
  - Negative 20 Degree Tilt
  - Right/Left Mouse Tray
- Warranty:
  - Mechanical – (3) Years
  - Electrical – (3) Years
  - Battery Performance – (30) Months



DCT-6LiA



DCT-6LiA-LCD

### Power System

Battery Type: 50ah Lithium Iron Phosphate

Inverter:

- USA – 150W, 120V 60Hz
- International – 230W, 230V 50/60Hz
- Certifications: UL60601-1 3rd Edition UL 1778 4th Edition



Power System Fuel Gauge

### Optional Accessories

- Slide-out Keyboard Tray
- Sani-Wipe, Scanner, Sharps Container, and Cup Holders
- Adjustable Laptop Screen Support
- Multiple VESA Mounts
- Several Camera Mounts
- Battery Power Packs
- Printer/Scanner/Equipment Shelves
- Storage Drawers with Lock Options, Wire Baskets
- Custom Paint Color