


Built from the ground up with medical-grade durability.

innovative-medical.com
114info@convergent.com
800.968.3135

Baby Bathinets



Bathinets bring the bath to the baby, creating a safer, convenient, and educational bathing experience within the patient room. Features include a washtub with intuitive temperature strip, heavy-duty soft-close drawers, smooth-rolling closed casters, and a protective bullnose bumper around the base. Available in a variety of colors, sizes, and configurations to support your workflow and room design.

Key benefits

Durable construction for the medical environment

Stands up to harsh daily use and repeated cleaning with all healthcare-grade solutions except 100% bleach (Antimicrobial only).

Unmatched customization

Our team works with you to recommend custom configurations that support your workflow, preferences, and space.

Backed by a 10-year warranty

Antimicrobial finishes are backed by a 10-year warranty for continued reliability (5-year on Laminate and Thermofoil).

We don't do "standard".

Our labor and delivery cart line is built to order, meaning the size, materials, and configuration are decided by you. See options below:

Overall size and dimensions	Color(s): See next page
Finish: Antimicrobial, Laminate, or Thermofoil	# of Drawers and/or cupboards
Flip-up work surface, built in handles	Locking on drawer(s) or cupboard(s)
Pass through drawers	Other customizations (just ask!)

White

Sanshade

Dolphin Gray

Nutmeg

Slate Gray

Charcoal

Indigo

Mocha Brown

Black

Finish and material shape daily cleaning practices and equipment lifespan.

01. Antimicrobial

Our most durable finish option, Antimicrobial is a recyclable polymer with inherent antimicrobial properties and color-through construction to help hide scratches. Compatible with all cleaning agents except 100% bleach, it meets Joint Commission standards, contributes to LEED points, and includes a **standard 10-year warranty**.



02. Laminate

Laminate offers a durable, cost-effective finish option with natural looking woodgrain finishes. Available in multiple grain patterns and colors, it includes a **standard 5-year warranty**.

03. Thermofoil

Thermofoil offers seamless edges and profiled design options for drawer and door fronts, though it is not as impact durable as laminate. Includes a **standard 5-year warranty**.

Kensington
Maple

Landmark
Wood
(Laminate Only)

Candlelight

Acajou
Mahogany

Chocolate
Pear

Walnut
Heights

Amber
Cherry

Laminate and Thermofoil Finishes

Trusted partner. Trusted solutions.

innovative-medical.com
114info@convergent.com
800.968.3135

convergent® | innovative
medical