

Built from the ground up with medical-grade durability.

innovative-medical.com
114info@convergent.com
800.968.3135

Bedside Cabinets



Medical bedside cabinets provide convenient bedside storage while helping simplify daily cleaning and disinfection. Smooth-rolling casters allow cabinets to move easily for access beneath and around the unit, supporting a cleaner patient environment. Available in a wide range of finishes, sizes, and configurations, these cabinets can be matched with bassinets, delivery carts, and other OB room furniture for a cohesive appearance throughout the space.

Key benefits

Durable construction for the medical environment

Stands up to harsh daily use and repeated cleaning with all healthcare-grade solutions except 100% bleach (Antimicrobial only).

Unmatched customization

Our team works with you to recommend custom configurations that support your workflow, preferences, and space.

Backed by a 10-year warranty

Antimicrobial finishes are backed by a 10-year warranty for continued reliability (5-year on Laminate and Thermofoil).

We don't do "standard".

Our labor and delivery cart line is built to order, meaning the size, materials, and configuration are decided by you. See options below:

Overall size and dimensions	Color(s): See next page
Finish: Antimicrobial, Laminate, or Thermofoil	# of Drawers and/or cupboards
Casters or no casters (mobile vs. stationary)	Locking on drawer(s) or cupboard(s)
Flip-up work surfaces, power strips, USB charging	Other customizations (just ask!)

White

Sanshade

Dolphin Gray

Nutmeg

Slate Gray

Charcoal

Indigo

Mocha Brown

Black

Finish and material shape daily cleaning practices and equipment lifespan.

01. Antimicrobial

Our most durable finish option, Antimicrobial is a recyclable polymer with inherent antimicrobial properties and color-through construction to help hide scratches. Compatible with all cleaning agents except 100% bleach, it meets Joint Commission standards, contributes to LEED points, and includes a **standard 10-year warranty**.



02. Laminate

Laminate offers a durable, cost-effective finish option with natural looking woodgrain finishes. Available in multiple grain patterns and colors, it includes a **standard 5-year warranty**.

03. Thermofoil

Thermofoil offers seamless edges and profiled design options for drawer and door fronts, though it is not as impact durable as laminate. Includes a **standard 5-year warranty**.

Kensington
Maple

Landmark
Wood
(Laminate Only)

Candlelight

Acajou
Mahogany

Chocolate
Pear

Walnut
Heights

Amber
Cherry

Laminate and Thermofoil Finishes

Trusted partner. Trusted solutions.

innovative-medical.com
114info@convergent.com
800.968.3135

convergent® | innovative
medical