

Labor & Delivery Carts:

Clinician and patient focused carts to support workflow, technology and aesthetics.

convergent[®] | innovative
medical

Convergent-Innovative Medical is experienced in designing carts that **perform reliably** against the demands of the 24/7/365 healthcare ecosystem. Whether standard, customized or fully custom, Convergent-IMS provides solutions for Labor & Delivery as well as other units in hospitals, clinics and extended care facilities.



Fetal Monitor



Bassinet



Delivery/Case



NICU



Overbed Table



Bedside



Custom

Features & Benefits



- ✓ **Workflow & Productivity**
 - **Display mounts:** Dual/single
 - **Flip-up side surface** for procedure kits or med prep
 - **Locks:** Keyed, punch code, electronic
 - **Optional accessories:**
 - Drawer dividers & lights
 - Power strips & cable wraps
 - Barcode scanner holder
 - Baskets
 - Keyboard tray: Slide-out, height adjustable



- ✓ **Infection Control**
 - **Tested** against common healthcare germicides
 - **Antimicrobial Finish** for infection control & durability
 - **Seamless Thermofoil finish** available on tops & drawer fronts for thorough cleaning
 - **Stainless steel handles** are standard
 - **SpillGuard tops** contain minor spills & small objects
 - **Premium casters** designed to minimize lint & hair accumulation



- ✓ **Ease of Use & Mobility**
 - **5" diameter premium closed casters** for ease of rolling
 - **Locking casters (2)**, optional directional (steering) caster available
 - **Detents in drawer & slideouts:** To keep 'closed' when cart is in motion
 - **Ergonomic-height** can be customized including handles on bassinet and keyboard tray on fetal monitor cart



- ✓ **Unsurpassed Durability**
 - Designed for the 24/7/365 demands of the acute healthcare environment
 - **Premium heavy-duty drawer slides** with full extension to reach back of drawers
 - **Bumper molding** around cart base helps protect the cart, walls and other equipment
 - **5-10 year Standard Warranties**