

Built from the ground up with medical-grade durability.

innovative-medical.com
114info@convergent.com
800.968.3135

Delivery Carts



Our Delivery Carts are built for durability, customization, and long-term performance in labor and delivery environments. Features include heavy-duty soft-close drawers, smooth-rolling closed casters that reduce hair and lint pickup, and a protective bullnose bumper around the cart base. Available in a wide range of finishes, sizes, and configurations, each cart is backed by a dependable warranty and tailored to your workflow needs.

Key benefits

Durable construction for the medical environment

Stands up to harsh daily use and repeated cleaning with all healthcare-grade solutions except 100% bleach (Antimicrobial only).

Unmatched customization

Our team works with you to recommend custom configurations that support your workflow, preferences, and space.


Backed by a 10-year warranty

Antimicrobial finishes are backed by a 10-year warranty for continued reliability (5-year on Laminate and Thermofoil).

We don't do "standard".

Our labor and delivery cart line is built to order, meaning the size, materials, and configuration are decided by you. See options below:

Overall size and dimensions	Color(s): See next page
Finish: Antimicrobial, Laminate, or Thermofoil	# of Drawers and/or cupboards
Flip-up side table(s) and wash basin cutouts	Locking on drawer(s) or cupboard(s)
Pass through drawers	Other customizations (just ask!)



Finish and material shape daily cleaning practices and equipment lifespan.

01. Antimicrobial

Our most durable finish option, Antimicrobial is a recyclable polymer with inherent antimicrobial properties and color-through construction to help hide scratches. Compatible with all cleaning agents except 100% bleach, it meets Joint Commission standards, contributes to LEED points, and includes a **standard 10-year warranty**.

02. Laminate

Laminate offers a durable, cost-effective finish option with natural looking woodgrain finishes. Available in multiple grain patterns and colors, it includes a **standard 5-year warranty**.

03. Thermofoil

Thermofoil offers seamless edges and profiled design options for drawer and door fronts, though it is not as impact durable as laminate. Includes a **standard 5-year warranty**.

04. Stainless Steel

Stainless steel offers a highly durable, easy-to-clean surface, though it may show fingerprints and scratches more easily. Includes a **standard 5-year warranty**.



Laminate and Thermofoil Finishes



Trusted partner. Trusted solutions.

innovative-medical.com
114info@convergent.com
800.968.3135

convergent® | innovative
medical