

Built from the ground up with medical-grade durability.

# Fetal Monitor Carts

innovative-medical.com  
114info@convergent.com  
800.968.3135



Designed to adapt to your workflow, our fetal monitor carts offer extensive customization options to support a wide range of labor and delivery environments. Features include single or dual display mounts, a slide-out keyboard tray, slide-out fetal monitor tray, heavy-duty soft-close drawers, smooth-rolling closed casters, and a protective bullnose bumper around the base. Available in a wide range of finishes, sizes, and configurations.

## Key benefits

### Durable construction for the medical environment

Stands up to harsh daily use and repeated cleaning with all healthcare-grade solutions except 100% bleach (Antimicrobial only).

### Unmatched customization

Our team works with you to recommend custom configurations that support your workflow, preferences, and space.


### Backed by a 10-year warranty

Antimicrobial finishes are backed by a 10-year warranty for continued reliability (5-year on Laminate and Thermofoil).

### We don't do "standard".

Our labor and delivery cart line is built to order, meaning the size, materials, and configuration are decided by you. See options below:

|                                                          |                                                       |
|----------------------------------------------------------|-------------------------------------------------------|
| Overall size and dimensions                              | Color(s): See next page                               |
| Finish: Antimicrobial, Laminate, or Thermofoil           | # of Drawers and/or cupboards, locking or non-locking |
| Slide-out or height adjustable keyboard tray             | Single, dual, horizontal or vertical display mounts   |
| Power strips, scanner mounts, CPU holders, baskets, etc. | Other customizations (just ask!)                      |



## Finish and material shape daily cleaning practices and equipment lifespan.

### 01. Antimicrobial

Our most durable finish option, Antimicrobial is a recyclable polymer with inherent antimicrobial properties and color-through construction to help hide scratches. Compatible with all cleaning agents except 100% bleach, it meets Joint Commission standards, contributes to LEED points, and includes a **standard 10-year warranty**.



### 02. Laminate

Laminate offers a durable, cost-effective finish option with natural looking woodgrain finishes. Available in multiple grain patterns and colors, it includes a **standard 5-year warranty**.

### 03. Thermofoil

Thermofoil offers seamless edges and profiled design options for drawer and door fronts, though it is not as impact durable as laminate. Includes a **standard 5-year warranty**.



Laminate and Thermofoil Finishes

Trusted partner. Trusted solutions.

innovative-medical.com  
114info@convergent.com  
800.968.3135

convergent® | innovative  
medical