

Built from the ground up with medical-grade durability.

innovative-medical.com
114info@convergent.com
800.968.3135

Height Adjustable Bassinets



Height adjustable bassinets provide added ergonomic support for caregivers and mothers while offering flexible configurations for different care environments. Features include 4 independent height adjustment levers, tub-tilt limiters to help prevent tub flip-out, ergonomic handle bars, tubs with flame-retardant mattresses, pass-through drawers, smooth-rolling closed casters, and a protective bullnose bumper around the base. Available in 6 finishes and several configurations.

Key benefits

Durable construction for the medical environment

Stands up to harsh daily use and repeated cleaning with many healthcare-grade solutions.

Ergonomic height adjustment

Four height adjustment handles located on each corner (used independently) for easy standing or bedside access.

Backed by a limited lifetime warranty

The pneumatic spring is backed by a 5-year warranty, while the remaining cart structure is covered by a limited lifetime warranty.

We don't do "standard".

Choose from multiple cart builds, finishes, and configurations to fit your workflow.

Options	Specifications
Color: See next page for choice of 6 colors	Surface dimensions: 36"W x 21.25"D
	Surface height: 30.5" - 40.25"
One or two pass-through drawers	Overall width: 37.5"
	Base dimensions: 35"W x 22"D
Pass-through work surface <u>OR</u> one or two flip-up work surfaces	Drawer dimensions: 22.5"W x 13.25"D x 6"H
	Mattress dimensions: 26"L x 12.5"W

Bright White

Eiderdown

American Natural

Blanca Elme

Aria

Walnut

Summer Flame

Finish and material shape daily cleaning practices and equipment lifespan.

Kydex® 3D Laminate

Our height adjustable bassinets use durable materials designed for repeated cleaning and daily clinical use.

Its KYDEX® 3D laminate is thicker and more impact resistant than standard thermofoil, helping the bassinet maintain its appearance and performance over time. The material also resists degradation from a wide range of hospital cleaners and disinfectants while meeting GREENGUARD® Gold certification standards.

Patented aluminum lift columns provide strength, stability, and smooth height adjustment for caregivers. The PETG bassinet tub offers increased impact resistance compared to traditional acrylic tubs.

For added confidence, the pneumatic lift system is backed by a **5-year warranty**, while the remaining cart structure is covered by a **limited lifetime warranty**.



Trusted partner. Trusted solutions.

innovative-medical.com
114info@convergent.com
800.968.3135

convergent® | innovative
medical