

Avalo® Series Medical Cart Line Anesthesia Cart

The Avalo Series Anesthesia Cart is designed with a synergy of function and technology to meet the unique needs of your anesthesia department. Keyless entry is standard on every Avalo Series Anesthesia Cart and the innovative guide system ensures drawers will not be left open. A flexible divider system, removable drawer trays, and a host of organizer accessories neatly arrange the medications and all related supply items.



Accessories & Options | Contact Us For More!



Bulk Storage Bins



Sharps Container




I.V. Pole



Divided Drawer Tray



Tilt Bin Organizer



Slide-Out Shelf

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Model Options



Standard



Intermediate



Compact

Anesthesia Cart Features

- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- Precision-welded steel frame with high-impact panels
- 5-year warranty
- Easy ordering, fast delivery

Specifications, availability & components are subject to change without notice.



Models and Components Specifications

Avalo® Series AC Anesthesia Cart



Model ACS



Model ACM

	Standard Height		Compact Height	
H	43"	110 cm	39.5"	100 cm
D	24"	61 cm	24"	61 cm
W	25.5"	65 cm	25.5"	65 cm

	Standard Height		Compact Height	
H	43"	110 cm	39.5"	100 cm
D	24"	61 cm	24"	61 cm
W	37.75"	96 cm	37.75"	96 cm



Main Drawers

	3" Main		6" Main		10" Main	
H	3"	7.6 cm	6.3"	16 cm	9.6"	24 cm
D	14.7"	37 cm	14.7"	37 cm	14.7"	37 cm
W	16.5"	42 cm	16.5"	42 cm	16.5"	42 cm




Utility Drawers

	3" Utility		6" Utility		10" Utility	
H	3"	7.6 cm	6.3"	16 cm	9.6"	24 cm
D	14.7"	37 cm	14.7"	37 cm	14.7"	37 cm
W	7.55"	19 cm	7.55"	19 cm	7.55"	19 cm



Cassettes

	2-Tier		3-Tier	
H	6.6"	16.8 cm	9.9"	25.1 cm
D	16.25"	41.3 cm	16.25"	41.3 cm
W	18"	45.7 cm	18"	45.7 cm



Bins

	8" Bin		5" Inch Bin	
H	2.75"	7 cm	2.75"	7 cm
D	12.5"	31.75 cm	12.5"	31.75 cm
W	8"	20.3 cm	5.25"	13.3 cm





Avalo® Series Medical Carts

Options & Accessories



Divided Drawer Tray



Drawer Dividers



Sharps/Glove Holder



Storage Container



Cardiac Board




I.V. Pole



Drawer Tag Holder



Oxygen Tank Holder



Accessory Bridge



Bulk Storage Bins



Tilt Bin Organizers



Tracking Wheel



Clear Top Mat



Storage Module



Internal Lock Box



Plug Outlet Strip
UL1363A Rated



Defibrillator Shelf



Waste Container



Writing Surface



Computer Arm