

Bassinets

Comfort and dependability for your smallest patients

convergent® | innovative medical

Standard. Twin. Compact. Our bassinets are customizable to meet the multiple demands from storing supplies for staff, mom and infant, and to provide peace of mind when infant is out of one's arms.



✓ Benefits

- Designed and built for the 24/7/365 demands of the healthcare environment
- Premium heavy-duty drawer slides with detents to prevent opening while cart is in motion
- Tested to withstand healthcare germicidal solutions
- Bullnose bumper around cart base
- Standard **5-Year Warranty** (Laminate, Thermofoil)
- Standard **10-Year Warranty** (Antimicrobial)

✓ Standard Dimensions/Configurations

- **34" W x 21" D x 37.5" H**
- 2 or 3 drawers
- Drawer over/under cupboard
- Custom Configurations

✓ Standard Features

- **FDA 510(k) registered:** 10054532, device D291909
- Dual Stainless Steel Handles are non-porous, infection control surfaces, extended for easy stride, ergonomically positioned
- Tub-tilt limiters prevent tub flip-out
- Dual reflex blocks to help "spitty" babies
- Premium, precision 5" diameter closed casters (front locking) for smooth, easy-rolling on carpet or hard flooring. Minimizes air and lint accumulation

✓ Available Options

- **Antimicrobial Finish** for infection control & durability (10y Warranty)
- **Thermofoil** tops and drawer fronts for a seamless finish (5y Warranty)
- Slide out tray on 100 lb. rated slides for equipment
- Locking drawers & doors
- Drawer dividers
- IV pole and mount
- Custom colors, designs, hardware, & accessories to accommodate your room design and workflow

Cleanability and Infection Control Design

✓ Choice of Finishes: Ask about antimicrobial color choices



Kensington Maple



Candlelight



Amber Cherry



Acajou Mahogany



Chocolate Pear Tree