

# Accessory Options for Altus Powered Workstations




**Small Basket**  
SKU: HB-PL  
Weight: 3 lbs.  
Holds bar code charging bases, zebra printers, etc. 12.28" W x 5"D x 4"H, column mounted, platinum (PL) finish



**Large Basket**  
SKU: HFH-PL  
Weight: 5 lbs.  
Holds various folders & files, 12" W x 4"D x 10 3/4"H, column mounted, platinum (PL) finish



**Rear Basket**  
SKU: RB-3-L  
Weight: 1 lbs. Stores blood pressure cuffs, IV bags, cordless barcode holder chargers, etc. 16.25"W x 5.24"D x 6"H, light grey (L) finish, weight limit not to exceed 5 lbs.



**Worksurface-Mounted Single Drawer, Locking**  
SKU: DSW  
Weight: 12 lbs. Drawer dimensions: 12.5"W x 2.29"H x 9 1/4"D worksurface mount, white finish, built in lock, includes 2 bins.




**Column-Mounted Single Drawer, Locking**  
SKU: DS1  
Weight: 11.5 lbs.  
Drawer dimensions: 9 1/4"D x 2.29"H x 12.5"W, white finish, built in lock, includes 2 bins: 4.35"W x 2.62"H x 12.372"D front label insert area fits business card size label each bin includes 2 adjustable bin dividers



**Column-Mounted Double Drawer, Locking**  
SKU: DS2  
Weight: 23 lbs.  
Drawer dimensions: 9 1/4"D x 2.29"H x 12.5"W, full outer dimensions: 10"D x 9 1/4"H x 14 1/2"W, white finish, two built in locks includes 4 bins: 4.35"W x 2.62"H x 12.372"D, includes 2 adjustable bin dividers



**Column-Mounted Single Drawer, Non-Locking**  
SKU: DS1-NL  
Weight: 11.5 lbs. Drawer dimensions: 9 1/4"D x 2.29"H x 12.5"W, white finish includes 2 bins: 4.35"W x 2.62"H x 12.372"D front label insert area fits business card size label each bin includes 2 adjustable bin dividers




**Worksurface-Mounted Single Drawer, Non-Locking**  
SKU: DSW-NL  
Weight: 12 lbs..  
Drawer dimensions: 12.5"W x 2.29"H x 9 1/4"D worksurface mount, white finish, includes 2 bins.



**Column-Mounted Double Drawer, Non-locking**  
 SKU: DS2-NL  
 Weight: 23 lbs.  
 Drawer dimensions: 9 1/4"D x 2.29"H x 12.5"W, full outer dimensions: 10"D x 9 1/4"H x 14 1/2"W, white finish includes 4 bins: 4.35"W x 2.62"H x 12.372"D, includes 2 adjustable bin dividers



**ClioAir Side Basket**  
 SKU: SB1  
 Weight: .5 lbs.  
 Stores blood pressure cuffs, IV bags, cordless barcode holder chargers, etc. 10" W X 4.75" D X 3.6" H, light grey (L) finish, weight limit not to exceed 5 lbs. Can be mounted on the left or right side of the workstation.



**ClioAir Rear Basket**  
 SKU: SB2  
 Weight: .5 lbs.  
 Stores blood pressure cuffs, IV bags, cordless barcode holder chargers, etc. 10" W X 4.75" D X 3.6" H, light grey (L) finish, weight limit not to exceed 5 lbs. Can be used on both the left and right side of the workstation.



**ClioAir Side Basket with Accessory Insert Mount**  
 SKU: SB3  
 Weight: .5 lbs.  
 The ClioAir side basket with an added accessory insert mount that acts as a false bottom with mounting holes. It mounts printers, bar code scanners, and more higher up for easy access. Can be used on both the left and right side of the workstation. 10" W X 4.75" D X 3.6" H, light grey (L) finish, weight limit not to exceed 5 lbs. Can be mounted on the left or right side of the workstation.



**Adjustable Bar Code Holder**  
 SKU: BCH-14-PL  
 Weight: 1.5 lbs.  
 Extends from 10"W to 12"W for right hand, left hand and top mount uses, VESA bracket monitor mount, platinum (PL) finish




**Monitor Arm Mounted Shelf**  
 SKU: BCHPS3  
 Weight: 1.5 lbs.  
 Includes mounting bracket and one full width shelf with slots for secure mounting. Each side 4.5"D x 9.25"W.



**Mouse House**  
 SKU: MH  
 Weight: .5 lbs  
 2.75"W x 4"H x 1.5"D Adhesive mount, white finish



**Column-Mount Brackets for Rear Mounting**  
 SKU: CMB  
 Weight: 1 lb  
 Attaches to the front of column to allow accessories to be mounted on the back side of the column. Front mount accessories can be used along with these brackets. 2 brackets (1 set), rear column mount, ionic white (IW) finish. If 2 rear mounted accessories are required, 2 sets of brackets need to be ordered along with the accessories. Accessories Sold Separately.



**Column-Mounted Sani-Wipe Bracket ONLY**  
 SKU: SW-1  
 Weight: 1 lb. Bracket allows attachment of sani-wipe mount & canister to Altus powered carts (designed to work with PDI 3-in-1 universal sani-bracket, #P010901) DOES NOT INCLUDE MOUNT TO HOLD CANISTER OR CANISTER  
 Bracket dimensions: 3.5" W x 2.8"H column mount, platinum (PL) finish



**Column-Mounted Sani-Wipe Bracket and Basket**  
 SKU: SW-2  
 Weight: 1 lb. Bracket allows for mounting of SaniWipe canister to the column of Altus workstation. It comes with a basket to secure the canister in place.




**Micro USB Powered Task Light**  
 SKU: TLKIT  
 Weight: .25 lbs.  
 Dimensions: 10" W x 3/8"H x 9/16"D, kit includes: 4, 5" zip ties and tie mounts, alcohol wipe, 10" L x .5"W velcro strip, 36" USB 2.0 male to micro B USB male white cable, anodized finish

For Additional Accessories, [Contact Convergent.](#)



**Multi-Purpose Adjustable Shelf**  
 SKU: MPS-6-PL  
 Weight: 6 lbs. Adjustable secure strap included 18"W x 7"-12"D x 2.6"H, column mount, platinum (PL) finish



**Cool Feet Laptop Risers**  
 SKU: CF  
 Weight: .5 lbs. Easily raise the height of the laptop in Altus laptop carts



**Electric Bonding**  
 SKU: EB10  
 Weight: .5 lbs. Electrically connects to any exposed metal that is not designed to carry electricity to protect from electric shock. Including: base, stretcher, weldment, handle, etc.

Accessory Compatibility Chart						
Accessories	Laptop Carts		LCD Carts			
		C1 / C1		C3/E3	L3 / L4	EL3 / EL4
BCH-14		o		x	x	x
BCHPS-3		o		x	x	x
CF		x		o	o	o
CMB		x		x	x	x
DSW / DSW-NL		x		x	o	o
EB10		x		x	x	x
MH		x		x	x	x
DS1 / DS2 / DS1-NL / DS2-NL		x		x	x	x
HB		x		x	x	x
HFH		x		x	x	x
MPS-6		x		x	x	x
SW-1		x		x	x	x
SW-2		x		x	x	x
RB-3-L		x		x	o	o
TLKIT		o		x	x	x
U222-004-R		x		x	x	x
SB1 / SB2 / SB3		o		o	x	x