

# Compact Bassinets

## Lower height for fall risk avoidance

convergent® | innovative medical

**Compact. Comfort. Fall-Risk.** Designed to work in conjunction with two prevalent practices: rooming-in and fall risk avoidance. The compact bassinet is lower, to afford mom an easier reach for the bassinet, whether she's had a C-section or difficult delivery.

Standard Height



Compact Height



### ✓ Benefits

- Designed and built for the 24/7/365 demands of the healthcare environment
- Premium heavy-duty drawer slides with detents to prevent opening while cart is in motion
- Tested to withstand healthcare germicidal solutions
- Bullnose bumper around cart base
- Standard **5-Year Warranty** (Laminate, Thermofoil)
- Standard **10-Year Warranty** (Antimicrobial)

### ✓ Standard Dimensions/Configurations

- **34" W x 21" D x 37.5" H**
- 1 or 2 drawers
- Custom Configurations

### ✓ Choice of Finishes: Ask about antimicrobial color choices

### ✓ Standard Features

- Tub-tilt limiters to prevent tub flip-out caused by sibling pulling up for better view
- Dual reflex blocks to help "spitty" babies
- Premium 5" diameter closed casters (2 locking/2 non-locking) for smooth, easy rolling, on carpet, or hard flooring. Minimizes hair and lint accumulation
- **Rooming-in Height** reduces strain on mom while reaching, reduces fall risk
- Dual handles of differing heights provide access to mom and clinician, are non-porous, infection control surfaces, extended for easy stride, ergonomically positioned

### ✓ Available Options

- **Antimicrobial Finish** for infection control & durability
- **Thermofoil Finish** tops and drawer fronts for a seamless finish
- Slide out tray on 100 lb. rated slides for equipment
- Locking drawers & doors
- Drawer dividers
- IV pole and mount
- Custom colors, designs, hardware, and accessories to accommodate your room design and workflow

Designed for Safety and Ergonomics



Kensington Maple



Candlelight



Amber Cherry



Acajou Mahogany



Chocolate Pear Tree