

Avalo® Series Medical Cart Line

Critical Care Cart

The Avalo Series Critical Care Cart is built to promote efficient workflow and task management in critical care environments. Designed for durability and flexibility, the Avalo Critical Care cart can be configured to meet your facility requirements.

The durable frame and spread caster position ensure stability, while making the cart easy to move and control. A choice of lock systems and full complement of accessory options permit customization to accommodate a variety of monitors, storage units and the security of critical care supplies.



Accessories & Options | Contact Us For More!



Bulk Storage Bins



Sharps Container




I.V. Pole



Waste Container



Tilt Bin Organizer



Slide-Out Shelf

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Model Options



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.

Critical Care Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty



Avalo® Series Medical Carts

Specifications

Standard
Height



Intermediate
Height



Compact
Height



Cart Specifications

	Surf. Height	Depth	Width
Standard Height:	43.0"	24"	31"
	110 cm	61 cm	79 cm
Intermediate Height:	39.5"	24"	31"
	100 cm	61 cm	79 cm
Compact Height:	36.0"	24"	31"
	91 cm	61 cm	79 cm

Drawer Specifications

Height	Depth	Width
3"	17.5"	21.5"
7 cm	45 cm	55 cm
6.25"	17.5"	21.5"
16 cm	45 cm	55 cm
10"	17.5"	21.5"
25 cm	45 cm	55 cm

Specifications, availability & components are subject to change without notice.

Avalo® Medical Carts are Certified



RoHS Restriction of
Hazardous Substances



ISO Quality
Management System
Standard



ISO Quality Management
System Standard



UL Meets UL Safety
Requirements



TUV Recognized for
Quality, Safety and
Sustainability

Avalo® Series Medical Carts

Options & Accessories



Divided Drawer Tray



Drawer Dividers



Sharps/Glove Holder



Storage Container



Cardiac Board




I.V. Pole



Drawer Tag Holder



Oxygen Tank Holder



Accessory Bridge



Bulk Storage Bins



Tilt Bin Organizers



Tracking Wheel



Clear Top Mat



Storage Module



Internal Lock Box



Plug Outlet Strip
UL1363A Rated



Defibrillator Shelf



Waste Container



Writing Surface



Computer Arm