

Avalo® Series Medical Cart Line

Emergency Cart

The Avalo Series Emergency Cart defines a new standard of organization and function, with a host of features that provide reliable performance and promote improved workflow in emergency departments.

The unique breakaway locking handle offers optimal security and the simple access essential to every emergency response. Stable maneuverability is built-in to every cart and a full line of accessories enhance user efficiency and support code team performance.



Quick-Access Lock

Avalo Emergency Cart's unique breakaway locking handle is standard on every model, offers the security you require and ensures quick access to code cart supplies. With a simple and swift motion, you can quickly break the tamper-evident seal and gain access to all drawers.

Emergency Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Quick-access lock offers security, yet provides quick access for code response
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty

Accessories & Options | Contact Us For More!



Sharps Container



Oxygen Tank Holder



I.V. Pole



Cardiac Board



Defibrillator Shelf



Slide-Out Shelf

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Model Options (also available in red)



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.



Avalo® Series Medical Carts

Specifications

Standard
Height



Intermediate
Height



Compact
Height



Cart Specifications

	Surf. Height	Depth	Width
Standard Height:	43.0"	24"	31"
	110 cm	61 cm	79 cm
Intermediate Height:	39.5"	24"	31"
	100 cm	61 cm	79 cm
Compact Height:	36.0"	24"	31"
	91 cm	61 cm	79 cm

Drawer Specifications

Height	Depth	Width
3"	17.5"	21.5"
7 cm	45 cm	55 cm
6.25"	17.5"	21.5"
16 cm	45 cm	55 cm
10"	17.5"	21.5"
25 cm	45 cm	55 cm

Specifications, availability & components are subject to change without notice.

Avalo® Medical Carts are Certified



RoHS Restriction of Hazardous Substances



ISO Quality Management System Standard



ISO Quality Management System Standard



UL Meets UL Safety Requirements



TUV Recognized for Quality, Safety and Sustainability



Avalo® Series Medical Carts

Options & Accessories



Divided Drawer Tray



Drawer Dividers



Sharps/Glove Holder



Storage Container



Cardiac Board




I.V. Pole



Drawer Tag Holder



Oxygen Tank Holder



Accessory Bridge



Bulk Storage Bins



Tilt Bin Organizers



Tracking Wheel



Clear Top Mat



Storage Module



Internal Lock Box



Plug Outlet Strip
UL1363A Rated



Defibrillator Shelf



Waste Container



Writing Surface



Computer Arm