

Fetal Monitor Carts

Carts that enhance your workflow

convergent[®] | innovative
medical

Ergonomic. Efficient. Enhancements. An unsung hero of carts, these support a host of equipment, including dual-monitors for charting, no matter the space available. Our carts also support your teams workflow, safety and ergonomic needs, and blend into your room setting unobtrusively.



✓ Benefits

- Designed and built for the 24/7/365 demands of the healthcare environment
- Premium heavy-duty drawer slides with detents to prevent opening while cart is in motion
- Tested to withstand healthcare germicidal solutions
- Bullnose bumper around cart base
- Standard **5-Year Warranty** (Laminate, Thermofoil)
- Standard **10-Year Warranty** (Antimicrobial)

✓ Choice of Finishes

- Ask about antimicrobial color choices

✓ Available Options

- **Antimicrobial Finish** for infection control & durability
- **Thermofoil Finish** tops and drawer fronts for infection control
- SpillGuard thermofoil top
- Optional accessories to customize your cart include barcode scanner mount, CPU holder, baskets, cord wrap, etc.
- Custom Designs, finishes hardware to accommodate your room design and workflow

✓ Display Mounting Options

- Single display, flush mount
- Dual stacked, flush mount
- Dual side-by-side with swivel or height-adjustable arms

✓ Workflow & Safety Options

- Slide out tray for keyboard
- Articulating/height adjustable arm for keyboard tray
- Slide-out tray for Fetal Monitor
- IV pole and mount
- Flip-up side table(s) add space for prep
- Locking drawers & doors
- Drawer dividers
- Multi-plug outlet strip - hospital grade

✓ Standard Dimensions/Configurations

- **25" W x 26" D x 40" H**
- 2 or 3 drawers available
- 2-drawer with hinged CPU compartment
- Drawer/Doors
- Custom Configurations



Kensington Maple



Candlelight



Amber Cherry



Acajou Mahogany



Chocolate Pear Tree