

# Isolation Cart

## Avalo® Series Medical Cart Line

The Avalo Series Isolation Cart provides a host of innovative features that promote long-term, reliable function and value. Large-volume seamless drawer trays are supported by a unique guide system that ensures consistent performance year after year.

Avalo's superior design offers durable service for long-term value and a solid feel uncommon in bent metal product alternatives.



### Isolation Cart Features

- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- Precision-welded steel frame with high-impact panels
- 5-year warranty
- Easy ordering, fast delivery


### Accessories & Options | Contact Us For More!



Keyless Locking\*



Best® Core Locks\*



Drawer Dividers



Clear Top Mat



Drawer Label Holder



Waste Container

### Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

### Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

### Model Options



Standard



Intermediate



Compact

*Specifications, availability & components are subject to change without notice.*



# Avalo® Series Medical Carts

## Specifications

Standard  
Height



Intermediate  
Height



Compact  
Height



### Cart Specifications

	Surf. Height	Depth	Width
<b>Standard Height:</b>	43.0"	24"	31"
	110 cm	61 cm	79 cm
<b>Intermediate Height:</b>	39.5"	24"	31"
	100 cm	61 cm	79 cm
<b>Compact Height:</b>	36.0"	24"	31"
	91 cm	61 cm	79 cm

### Drawer Specifications

Height	Depth	Width
3"	17.5"	21.5"
7 cm	45 cm	55 cm
6.25"	17.5"	21.5"
16 cm	45 cm	55 cm
10"	17.5"	21.5"
25 cm	45 cm	55 cm

Specifications, availability & components are subject to change without notice.

### Avalo® Medical Carts are Certified



**RoHS** Restriction of Hazardous Substances



**ISO** Quality Management System Standard



**ISO** Quality Management System Standard



**UL** Meets UL Safety Requirements



**TUV** Recognized for Quality, Safety and Sustainability



# Avalo® Series Medical Carts

## Options & Accessories



Divided Drawer Tray



Drawer Dividers



Sharps/Glove Holder



Storage Container



Cardiac Board




I.V. Pole



Drawer Tag Holder



Oxygen Tank Holder



Accessory Bridge



Bulk Storage Bins



Tilt Bin Organizers



Tracking Wheel



Clear Top Mat



Storage Module



Internal Lock Box



Plug Outlet Strip  
*UL1363A Rated*



Defibrillator Shelf



Waste Container



Writing Surface



Computer Arm