

Avalo® Series Medical Cart Line

I.V. Therapy Cart

Organizing I.V. preparation supplies is easier with the thoughtful design of the Avalo Series I.V. Therapy Cart. Optional drawer divider systems and storage bins offer the flexibility to organize and store I.V. supplies.

A large, sturdy work surface provides ample room to prepare I.V. solutions. A slide-out work surface extends to the left or right side when additional workspace is needed. A stable, adjustable I.V. pole can be positioned on the cart to meet your workflow needs.



Accessories & Options | Contact Us For More!



Drawer Dividers



Sharps Container




I.V. Pole



Waste Container



Side Organizer



Slide-Out Shelf

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Model Options



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.

I.V. Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty



Avalo® Series Medical Carts

Specifications

Standard
Height



Intermediate
Height



Compact
Height



Cart Specifications

	Surf. Height	Depth	Width
Standard Height:	43.0"	24"	31"
	110 cm	61 cm	79 cm
Intermediate Height:	39.5"	24"	31"
	100 cm	61 cm	79 cm
Compact Height:	36.0"	24"	31"
	91 cm	61 cm	79 cm

Drawer Specifications

Height	Depth	Width
3"	17.5"	21.5"
7 cm	45 cm	55 cm
6.25"	17.5"	21.5"
16 cm	45 cm	55 cm
10"	17.5"	21.5"
25 cm	45 cm	55 cm

Specifications, availability & components are subject to change without notice.

Avalo® Medical Carts are Certified



RoHS Restriction of Hazardous Substances



ISO Quality Management System Standard



ISO Quality Management System Standard



UL Meets UL Safety Requirements



TUV Recognized for Quality, Safety and Sustainability



Avalo® Series Medical Carts

Options & Accessories



Divided Drawer Tray



Drawer Dividers



Sharps/Glove Holder



Storage Container



Cardiac Board




I.V. Pole



Drawer Tag Holder



Oxygen Tank Holder



Accessory Bridge



Bulk Storage Bins



Tilt Bin Organizers



Tracking Wheel



Clear Top Mat



Storage Module



Internal Lock Box



Plug Outlet Strip
UL1363A Rated



Defibrillator Shelf



Waste Container



Writing Surface



Computer Arm