

# Scale Carts

## A mobile weighing cart for rooming-in

convergent<sup>®</sup> | innovative  
medical

**Mobile. Safety. Customizable.** Scale carts are used when rooming-in is strongly practiced. In some cases, the nursery has been repurposed, requiring in-room weighing of newborns. Having a mobile scale cart that looks more like the room's furniture is more appealing to family members and staff, than a cold, clinical stainless steel traditional scale cart.



### ✓ Benefits

- Designed and built for the 24/7/365 demands of the healthcare environment
- Premium heavy-duty drawer slides with detents to prevent opening while cart is in motion
- Tested to withstand healthcare germicidal solutions
- Bullnose bumper around cart base
- Standard **5-Year Warranty** (Laminate, Thermofoil)
- Standard **10-Year Warranty** (Antimicrobial)

### ✓ Choice of Finishes:

- Ask about antimicrobial color choices

### ✓ Standard Features

- Premium, precision 5" diameter closed casters for easy-rolling, even on carpet
- 2-locking and 2-non-locking
- Long lasting, minimizes hair and lint accumulation

### ✓ Standard Dimensions/Configurations

- **27" W x 22.25" D x 37.5" H**
- 2 or 3 drawer
- Drawer over/under cupboard
- Custom Configurations

### ✓ Workflow & Safety Options

- Single or dual flip-up side surfaces
- Galley rail on four-sides
- Locking drawers & doors
- Drawer dividers

### ✓ Available Options

- **Antimicrobial Finish** for infection control & durability
- **Thermofoil** tops and drawer fronts for a seamless finish
- Custom designs, finishes, hardware to accommodate your room design and workflow



Kensington Maple



Candlelight



Amber Cherry



Acajou Mahogany



Chocolate Pear Tree