

# Twin Bassinets

## Bassinets that deliver twice the joy

convergent<sup>®</sup> | innovative  
medical

**Customizable. Quiet. Safe.** Bassinets are customizable to meet the demands of twin births, by storing twice the supplies, keeping infants safe and secure, and providing peace of mind to mom and clinicians.



### ✓ Benefits

- Designed and built for the 24/7/365 demands of the healthcare environment
- Premium heavy-duty drawer slides with detents to prevent opening while cart is in motion
- Tested to withstand healthcare germicidal solutions
- Bullnose bumper around cart base
- Standard **5-Year Warranty** (Laminate, Thermofoil)
- Standard **10-Year Warranty** (Antimicrobial)

### ✓ Choice of Finishes

- Ask about antimicrobial color choices

### ✓ Standard Features

- Ergonomic, dual stainless steel transport handles are extended for easier stride
- Tub-tilt limiters to prevent tub flip-out caused by sibling pulling up for better view
- Dual reflex blocks to help "spitty" babies
- Premium 5" diameter closed casters (front locking) for smooth, easy rolling, on carpet, or hard flooring. Minimizes hair and lint accumulation

### ✓ Standard Dimensions/Configurations

- **42.025" W x 32.25" D x 37" H**
- 2 or 3 drawers
- Drawers over/under cupboard
- Custom Configurations

### ✓ Available Options

- **Antimicrobial Finish** for infection control & durability
- **Thermofoil Finish** tops and drawer fronts for a seamless finish
- **Workflow** and safety options:
  - Slide out tray on 100 lb. rated slides for equipment
  - Optional accessories to customize your cart
  - Locking drawers & doors
  - Drawer dividers
  - IV pole and mount
- Custom designs, finishes hardware to accommodate your room design and workflow



Kensington Maple



Candlelight



Amber Cherry



Acajou Mahogany



Chocolate Pear Tree