

Antimicrobial Carts:

Supporting the mixed demands of equipment and aesthetics

convergent[®] | innovative
medical

The labor and delivery environment is expected to “look like home” as well as manage the ever-growing world of technology. Offering **solutions that perform reliably** against the demands of the 24/7/365 healthcare ecosystem and designed to **enhance your workflow** to minimize or eliminate workarounds is key to offering a best-in-class solution.



Features & Benefits

- ✓ Exclusive **antimicrobial material** protects against surface stains and odor causing bacteria, algae and fungi
- ✓ **Scratch resistant**, color-through enhances aesthetics and infection control
- ✓ **Easy to clean**
- ✓ **Tested** to withstand healthcare germicidal solutions
- ✓ **Extended 10-year warranty**

Infection Control




- **Germicidal solutions:** Withstands all common wipes & solutions
- **Solid color-through finish:** Scratch resistant, will not delaminate, rot, chip or swell
- **Spill guard tops:** Available to contain minor spills
- **Premium dual casters:** Designed to minimize lint & hair accumulation, roll smoothly

Ease of Use & Mobility



- **5" diameter premium closed casters:** For ease of rolling, even on carpet
- **Locking casters:** (2) standard; optional directional (steering) caster available
- **Detents In drawer & slide-outs:** Keeps drawer ‘closed’ when cart is in motion
- **Ergonomic-height:** Handles on bassinet and keyboard tray on fetal monitor cart for majority population

Unsurpassed Durability



- **Designed and built for 24/7/365** demands of the acute healthcare environment
- **Premium heavy-duty drawer slides** with full extension to reach back of drawers
- **Bumper molding** around cart base helps protect the cart, walls and other equipment

Go Green! With Antimicrobial Finishes

convergent® | innovative
medical



Green qualifying polymer sheets are produced using processes that are non-polluting & conserve energy and natural resources. It is **economically and environmentally sound, safe and recyclable**. No harmful dusts are created during cutting and fabricating, no carcinogens or toxins. Also no gasses are emitted if burned, reducing the impact on our overall environment.

FDA-approved polymers are highly durable, ensuring quality and consistency in the final product. **Antimicrobial Polymer** provides a super tough surface for the many applications of the healthcare environment.

- Non-polluting
- No acids or harsh chemical for cleaning
- Completely recyclable
- Made from FDA-approved polymers
- Will not warp, rot or delaminate
- Antimicrobial against bacteria, algae & fungi
- Repels moisture
- Low maintenance
- Scratch resistant
- High impact
- Odor and stain resistant
- Chemical resistant
- Environmentally green
- Qualified for LEED Points
- Available in multiple colors

LEED Rating Points:

✓ Existing Buildings, September 2008

1 point: MR Credit 2.2 Sustainable Purchasing: Durable Goods

✓ New Construction, Version 2.2

2.2 points: MR credit 4.2 Recycled Content

✓ Commercial Interiors, Version 2.0

2 points: MR credit 4.2 Recycled Content

1 Point: EQ credit 4.5, Low-Emitting Materials

✓ Homes, January 2008

.5 Point: MR credit 2.2 Environmentally Preferable Products

✓ Schools, First Edition, November 2007

2 points: MR Credit 4.2 Recycled Content

✓ Green Guide for Health Care, Version 2.2

1 point: MR Credit 5.2 Furniture & Medical Furnishings: Materials

1 point: MR Credit 5.3 Furniture & Medical Furnishings: Manufacturing, Transportation and Recycling

Call or Email your Convergent-IMS rep to see how we can help your hospital Go Green!