

Supply Carts

Uniquely customized for any purpose



Functional. Customizable. Personal. Supply cart needs vary from one unit to another, as much as from one hospital to another. Maximize the utility by configuring a cart that meets your specific requirements, especially when it comes to blending it into the room.



✓ Benefits

- Designed and built for the 24/7/365 demands of the healthcare environment
- Premium heavy-duty drawer slides with detents to prevent opening while cart is in motion
- Tested to withstand healthcare germicidal solutions
- Bullnose bumper around cart base
- Standard **5-Year Warranty** (Laminate, Thermofoil)
- Standard **10-Year Warranty** (Antimicrobial)

✓ Available Options

- **Antimicrobial Finish** for infection control & durability
- **Thermofoil Finish** tops and drawer fronts for seamless design
- SpillGuard thermofoil top
- Optional accessories to customize your cart include barcode scanner mount, CPU holder, baskets, glove box, cord wrap, etc.
- Custom designs, finishes hardware to accommodate your room design and workflow

✓ Workflow & Safety Options

- Flip-up side table(s) add space for prep
- Cabinet vent if using as soiled linen cart
- Locking drawers & doors
- Drawer dividers
- Multi-plug outlet strip - hospital grade


✓ Display Mounting Options

- Single display, flush mount
- Dual stacked, flush mount
- Dual side-by-side with swivel or height-adjustable arms

✓ Standard Dimensions/Configurations

- **27" W x 26" D x 40" H**
- 4 or 5 drawers
- Custom Configurations

✓ Choice of Antimicrobial Finishes: Ask about laminate color choices



White

Dolphin Gray

Sansshade

Nutmeg

Mocha Brown

Slate Gray

Charcoal

Indigo

Black