

Built from the ground up with medical-grade durability.

innovative-medical.com
114info@convergent.com
800.968.3135

Overbed Tables

By MedViron



Overbed tables by MedViron provide a stable, adjustable work surface designed for patient rooms in clinical environments. Height adjustable from 33.5" to 45.5" with a 100 lb. weight capacity, these tables support meals, personal items, documentation, and bedside care activities. Anodized aluminum columns deliver strength and smooth height adjustment. Backed by a limited lifetime warranty with a 5-year warranty on the gas spring. Available in a wide range of colors, sizes, and configurations.

Key benefits

Durable construction for the medical environment

Stands up to harsh daily use and repeated cleaning with many healthcare-grade solutions.

Ergonomic height adjustment

Flexible height adjustment supports both patient needs and caregiver ergonomics.

Backed by a limited lifetime warranty

The pneumatic spring is backed by a 5-year warranty, while the remaining cart structure is covered by a limited lifetime warranty.

We don't do "standard".

Our labor and delivery cart line is built to order, meaning the size, materials, and configuration are decided by you. See options below:

Overall size and base profile	Color(s): See next page
Tablet holder or vanity mirror	Spillguard edge
Flip-up work surface	Cupholder cutouts
Hanging single or dual side shelves	Other customizations (just ask!)

Bright White

Eiderdown

American Natural

Blanca Elme

Aria

Walnut

Summer Flame

Finish and material shape daily cleaning practices and equipment lifespan.

Kydex® 3D Laminate

Our height adjustable overbed tables by MedViron use durable materials designed for repeated cleaning and daily clinical use.

Its KYDEX® 3D laminate is thicker and more impact resistant than standard thermofoil, helping the overbed table maintain its appearance and performance over time. The material also resists degradation from a wide range of hospital cleaners and disinfectants while meeting GREENGUARD® Gold certification standards.

Patented aluminum lift columns provide strength, stability, and smooth height adjustment for caregivers.

For added confidence, the pneumatic lift system is backed by a **5-year warranty**, while the remaining cart structure is covered by a **limited lifetime warranty**.



Trusted partner. Trusted solutions.

innovative-medical.com
114info@convergent.com
800.968.3135

convergent® | innovative
medical